

- 3 cuff-sizes available (standard, XL, XS)
- softmatic measurement with optimal inflation and fast deflation
- additional average value indicator for a better blood pressure control
- memory of 30 measurements for 14-day-profile
- irregular heart beat indication
- clock/acoustic reminder
- AC adapter available as optional accessory
- extra-large 3-value-display
- with carrying case, batteries and blood pressure record card
- distinguished with the Quality seal by The German Hypertension Society



Attachment for additional cuff sizes

All instruments come with a standard cuff. The correct cuff size is important for reliable results. Therefore, boso offers additional cuff sizes for upper arm instruments if required:
 XS = 16 – 21 cm | Standard = 22 – 32 cm
 Universal = 22 – 42 cm | XL = 33 – 48 cm

The comfortable blood-pressure-instrument with exceptional features. Fully automatic blood pressure measurement at the wrist.

- irregular heart beat indication
- evaluation of blood-pressure in accordance with WHO
- additional average value indicator for a better blood-pressure control
- memory for 30 measurements for 14-day-profile
- quick and gentle measurement during inflation
- extra-large 3-value-display
- with carrying case, batteries and blood pressure record card
- distinguished with the Quality seal by The German Hypertension Society



NL NEOLAB SA

Mezzi ausiliari per l'indipendenza a domicilio
 Forniture ospedaliere e per case anziani
 CH-6883 Novazzano
 Tel. 091 683 03 51
www.neolab.ch – info@neolab.ch