

SEDIA ROTELLE SURACE 200



Doppia crociera
 Braccioli desk ribaltabili
 Pedane girevoli ed estraibili
 Freni

DATI TECNICI

 cm Ø 30	 cm Ø 20	 L 43	 L 53	 L 64 - 85	 h 51
200 Large:		L 48	L 58		
200 Xlarge:		L 53	L 63		
200 XXlarge:		L 60	L 70		
 Standard	 Standard	 Standard	 Optional	 Standard	 Portata 130 kg
 Kg 14,5	=	 Kg 13	+	 Kg 1,5	200 Large: 140 kg 200 Xlarge: 150 kg 200 XXlarge: 160 kg