

# seca 760

Mechanical personal scales in modern colours

new



- Robust construction.
- Large, easy-to-read dials.
- Low profile and wide platform.
- Precise and durable.

## seca 760:

A classic in several colors.

This seca classic has greatly influenced the image of scales. With an absolutely precise dial, synthetic covering used in medical applications, and steel plate housing from the good old days. The modern age has also left its mark on these scales, which are available in many different colors and models. For one thing a classic should never be boring.

In several colors to suit every taste.

These scales are available in burgundy red, midnight blue and ecru. And because there's no disputing taste, this classic model also comes in two additional versions – with white or black housing. All of them have

in common the chrome-plated safety ring, a platform of skin-friendly grooved natural rubber, the classic scale dial and robust, precise weighing mechanics.



Traditional weighing mechanics in classic colors.

## seca 760

### Technical Data

- Capacity: 150 kg
- Graduation: 1.000 g
- Dimensions (WxHxD):  
300 x 105 x 460 mm
- Weight: 4 kg

